



## Load type and process/autoclave selection

*PDA Europe*

*Basics and Recent Developments for Autoclaves and  
Lyophilizers - Tel Aviv/Israel, 13 November 2017*

*Sara Ferretti  
Fedegari Group*

## Solid loads

### Glassware

### Metallic items

- Surgical instruments
- machine parts

### Porous loads

- filters
- textiles
- stoppers in bags (or not),



## Liquid loads

### Sealed containers (LVP, SVP)

- standard closed containers (bottles, ampoules, vials)
- variable volume containers (syringes)



### Not-sealed containers

- flask with culture medium
- carboy (with filter on vent)




Saturated steam autoclave

Superheated water autoclave

Steam-air mixture autoclave

```
graph LR; A([Counterpressure (moist heat) autoclaves]) --> B[Superheated water autoclave]; A --> C[Steam-air mixture autoclave];
```

Counterpressure  
(moist heat) autoclaves



Direct contact



*STEAM-MO  
CONTACT*



Indirect contact

Heating &  
Sterilizing



*STEAM  
FUNCTION*




Heating

What is the  
sterilizing agent?

# SATURATED STEAM AUTOCLAVE

The temperature and pressure of saturated steam have a one-to-one correspondence



Saturated steam curve:  
water vapour in equilibrium with liquid water at the same T

If you choose the temperature, the pressure inside the chamber is automatically determined!

Temperature and pressure inside the sterilizer chamber

T (° C)		P (bar abs)
110	→	1.5
121	→	2.05
135	→	3.1

Taking into account steam sterilization requirements...


Steam-MO contact




Air must be removed  
steam without non-condensable gases

They **stratify** at the bottom of the chamber because they are more dense than steam; they **limit the heat exchange** between steam and product

The difference between the two cycles is the method used to remove air (steam injection or steam-vacuum pulses)



Metal items, empty glassware..



Porous solids (where air removal is critical): filters, textiles, stoppers in bag (or not), hollow materials...

- ✓ **Steam** is fed rapidly into the chamber until the sterilization temperature is reached
- ✓ **Temperature equilibration/penetration time:** delay between temperature of the product (product probes) and temperature of the chamber (monitoring probes)
- ✓ The **condensate is continuously removed** by a flow of **dynamic steam**: the vacuum pump always extracts condensate through a small valve
  - ✓ **Fresh vapor** continuously replaces the removed steam
  - ✓ Excellent **stability** and **uniformity of temperature** inside the chamber

The selection depends on the load type and on the final required results (e.g. wet or dry product)

**Drying and “natural” cooling** by final vacuum (solids or very small sealed containers with liquids, e.g. ampoules)

**Indirect cooling** by cold water circulation in the jacket and/or in internal heat-exchangers (plates); with air counterpressure (liquids in non-sealed containers)

**Direct cooling** by water spray onto the load:

- with air counterpressure (liquids in sealed small containers, SVP)
- without air counterpressure (only for liquids in very small sealed containers)



A typical example of a porous solid load to be treated  
in a saturated steam autoclave



- ✓ Load: *Rubber Stoppers in Tyvek/Plastic bag*
- ✓ Autoclave type: *Saturated Steam Sterilizer*
- ✓ Requirements: *Residual Humidity  $\leq 0.1\%$*

## **Residual humidity $\leq 0,1 \%$**

- ✓ Strongly needed for freeze dried products
- ✓ no risk of microbial growth
- ✓ to preserve machinability
- ✓ to preserve packaging barrier capability

## **Humidity in rubber stoppers: potential sources**

- ✓ Stopper formulation/matrix
- ✓ Storage conditions (ex. in the sterile area at the filling facility)
- ✓ Sterilization treatment and further processing steps

Before the sterilization phase...


- ✓ Hot air is inserted into the chamber to heat the load and, therefore, to reduce the creation of condensate
- ✓ The removal of the air is carried-out with a modulated depressurization (to maintain the integrity of the bags) followed by modulated steam-vacuum pulses; the sequence (pulse) is repeated several times




Air removal by **steam injection**:  
Metal items, glassware

Air removal by **vacuum/steam pulses**:  
porous solids (where air removal is critical)


**STOPPERS  
IN BAGS**



Cycle with modulated vacuum/steam pulses to not damage the load (e.g. for filters, membranes, stoppers in bags...)




## Stoppers in bags: typical cycle



After the sterilization phase...

- ✓ The vacuum causes the rapid evaporation of the condensate deposited on the material; in order to evaporate, this condensate requires vaporization calories which it draws from the material, which therefore cools
- ✓ During these phases the jacket is full of steam, and the heat which is emitted contributes to evaporate the condensate
- ✓ Air is injected into the chamber to help the thermal exchange (the transmission of heat from the jacket is poor because of vacuum)

## Stoppers in bags: typical cycle



**Result:**  
**RH  $\approx$  0,015%**

- ✓ Material design & packaging system
- ✓ Item orientation/arrangement
- ✓ Load initial temperature  
(& sterilization temperature)
- ✓ Drying by vacuum improved by:
  - ✓ Auxiliary heating equipment
  - ✓ Vacuum/ (hot) air pulses
  - ✓ Forced circulation of hot air (ex. fan)





✓ Our results after cycle optimization...

Drying phase	Cycle time	Residual Humidity [%]
Jacket (steam) + hot air	298 min	0,228
Jacket (steam) + hot air + FAN	291 min	0,015

## A FOF autoclave could be the right choice

Compared to LVP, they are more resistant to pressure differences: the resistance decreases as the diameter increases



What are the other options to perform the cooling ?



**Indirect cooling**

by cold water circulation in  
the jacket and/or in plates  
with air counterpressure



**Direct cooling**

by water spray with  
air counterpressure

The choice depends on customer needs (i.e. cycle time, final unloading temperature, product unloaded wet or dry) !

## SOLIDS

### Hard/Porous loads

✓ Glassware, plastic tools (empty)



✓ Metallic items (machine components, surgical instruments, tools)



✓ Filters



✓ Textiles




✓ Stoppers in bag (or not)

✓ Wrapped items (steam permeable wrapping)



## LIQUIDS

✓ Culture media  
(not sealed containers)



✓ Glass ampoules



They resist to the  
overpressure!

✓ Glass vials

# COUNTERPRESSURE AUTOCLAVES

## PHARMA INDUSTRIES



## FOOD INDUSTRIES



The total pressure (P) generated inside the sealed container at the temperature T (ex. 121 °C) is equal to:

$$P = P_{v(T)} + P_{a(T)}$$

Pressure of the water  
vapour


Pressure of the air

1) air initially present in the head space;  
2) dissolved gases that come out of the solution;  
3) reduction of the head space due to the thermal expansion of the liquid

\*Sealed= hermetically closed

Bottle partially filled with water solution at 20°C and 1 bar abs  
subjected to a saturated steam sterilization at 121°C

Total pressure in the chamber:  
Partial steam pressure= 2,05 bar



Partial steam pressure = 2,05 bar  
Partial air pressure= 1,34 bar

---

Total pressure in the bottle  $\cong$  3,4 bar

$$\Delta P \cong 1,4 \text{ bar} \rightarrow 1,4 \text{ kg/cm}^2$$

The total pressure inside the chamber is automatically controlled and adjusted according to:

- ✓ **Temperature of the solution**

- ✓ **Container features**  
(ex. rigid or deformable material)


✓ Steam-air mixture  
autoclaves

- Should the load be unloaded dry?
- Should the load temperature transition be fast?



✓ Superheated  
water  
autoclaves






Cooling under pressure is always used but...

- ✓ Superheated water autoclave: **DIRECT** cooling by water spray
  - Product is unloaded **wet**
  - The cooling phase is **faster** (heat exchange occurs through a liquid: more efficient)
  
- ✓ Steam-Air mixture autoclave: **INDIRECT** cooling by forced air circulation (fan) + cold water circulation in the plates and jacket (if present)
  - Product can be unloaded **dry**
  - The cooling phase is **longer** (heat exchange occurs through a gas: less efficient)

The main issues to consider while sterilizing PFSs with a moist heat process are:




How to avoid  
**plunger expulsion**  
during treatment



The head space and the  
space behind the plunger  
is small



## Container sealed with a plunger: variable volume container



$Q$  = amount of liquid

$V_1$  = head space volume for container 1

$V_2$  = head space volume for container 2

**PRESSURE on PLUNGER 1 > PRESSURE on PLUNGER 2**

An aqueous solution increases its volume about 6%  
when heated from ambient temperature to 121°C  
Therefore...



Thermal expansion of the water becomes important if  
the head space is lower than 10-15% of the volume of  
the container



There is **no practical mean to prevent the thermal expansion of liquids** during sterilization

The pressure required to reduce of 6% the volume of a liquid like water would be very large: thousands of bars !!!



**FAQ**

- ✓ Is the material resistant to high temperature?
- ✓ What is the melting point of the material?

The answers can tell us if steam sterilization is  
the right choice and how to develop the  
sterilization process

What happen if...


- ✓ Temperature is too high?
- ✓ The applied counterpressure is not well adjusted?




**Container  
deformation**



Temperature probe  
inside the sample



## *Rotating basket for liquids. Why??*





Emulsions,  
Suspensions



Heat sensitive  
products  
(sometimes)



Dense or non  
homogeneous  
mixtures

### Superheated water autoclaves

- ✓ Easy control of heating and cooling rate
- ✓ Short process duration
- ✓ No consumption of clean steam (used only for filter sterilization)
- Product is unloaded wet
- Higher water consumption (for initial filling)
- Blushing phenomenon (i.e. whitening of the PVC due to water absorption)

### Steam-air mixture autoclaves

- Indirect and 'difficult' control of heating and cooling rate
- Longer process duration (mainly because of indirect cooling)
- Consumption of clean steam
- ✓ Product could be easily unloaded dry
- ✓ No PW/UPW/WFI water consumption
- ✓ Blushing phenomenon very<sup>49</sup> rare

If your load is **compatible with moist heat sterilization** conditions (ex. temperature, pressure and humidity) you must use this method to sterilize it

Is your **load resistant to overpressure?**

Yes → Saturated steam autoclaves

No → Counterpressure steam autoclaves

- ✓ Is the product **solid** or **liquid**?
- ✓ What are the characteristics of the **container**?
  - Plastic or glass material
  - Sealed container (ex. «standard» sealing or variable volume container) or open container
  - Large or small volume
- ✓ Should the product be unloaded **dry** or **wet**?
- ✓ Does the product need to be **rotated**?

## LIQUIDS

Open container

Sealed container

**Saturated steam  
autoclaves**

SVP  
( $\leq 100$  mL)


LVP  
( $> 100$  mL)

Non deformable  
container

Deformable  
container

**Counterpressure  
autoclaves**

SVP= Small Volume Parentals  
LVP= Large Volume Parentals





Thank you

saf@fedegari.com  
training@fedegari.com