



# ControLyo™ Technology and TDLAS in Commercial Manufacturing as a QBD tool supporting scale up

*Application of scalable tools to aid process development in scale up, and batch recovery by scale down*

*PDA Europe*

*Event*

*City/Country, Day Month Year*

*Ian Whitehall*

*CMO*

*SP Scientific*


# ICH Q11 - Development and Manufacture of Drug Substances

Summary Statement of ICH Q11 guidance:

- Identifying potential CQAs associated with the drug substance so that those characteristics having an impact on drug product quality can be studied and controlled
- Defining an appropriate manufacturing process
- Defining a control strategy to ensure process performance and drug substance quality

## Importance of PAT

- Enhance understanding of critical product attributes which can impact over final quality
- Characterize product temperature profile and product resistance during development and scale up
- Rationalize information in product life cycle management and quality decision
- Regulatory expectation - *PAT – A Framework for Innovative Pharmaceutical Development, Manufacturing and Quality Assurance (FDA Guidance for Industry)*



# ControlLy<sup>TM</sup> Technology Manufacturing Adoption


- 40 SP LyoStar3 R&D freeze dryers with ControlLy<sup>TM</sup> technology in use world wide

## Commercial Units:

- Evaluation of production batch for stability study in human injectable products
- Use in commercial manufacturing for animal health products


# Scanning Electron Microscope (SEM) Image of Uncontrolled and Controlled Freezing

**Uncontrolled Freezing**



**Using 1°C/min shelf cooling rate**

**ControlLyo™ Freezing -3°C**



**Using ControlLyo™ at -3°C Shelf SP**

**SEM images of sucrose, 75 mg/mL**

SEM Image

# **ControlLy<sup>™</sup> in a Manufacturing Environment**

- Scalability in freeze dryers of any size
- Increases product consistency and uniformity
- Less vial damage
- Reduces cycle times and improves product yield
- Technology Differentiator – added capabilities
- Robust, non-invasive, and easily implemented/maintained
- Conforms to regulatory expectation

## Recent Work Using **ControlLyo™**


- “Application of **controlled nucleation** during lyophilization to improve cake appearance and product quality” – SP Webinar by Dr. Stuart Wang, (formerly of Biogen IDEC)
- “Impact of **controlled ice nucleation** on process performance and quality attributes of a lyophilized monoclonal antibody”, FDA, Awotwe-Otoo, D., Agarabi, C., Read, E., Lute, S., & Borson, K. (2013), *International Journal of Pharmaceutics*, 450(1-2), 70–78.




# PAT: TDLAS Mass Flow Measurements

Near-IR ( $\sim 1.4 \mu\text{m}$ ) water vapor absorption measurements to determine:

- 1) Gas temperature (K)
  - 2) Water vapor concentration [molecules/ $\text{cm}^3$ ]
  - 3) Gas flow velocity [m/s]
- Calculate the water vapor flow rate,  $dm/dt$  [grams/s]
- Integrate the water removal rate to predict the mass balance



Absorption lineshapes from two line-of-sight measurements across the spool connecting the lyophilizer chamber and condenser




# TDLAS Measurement Applications

- Determination of primary and secondary drying endpoints
- Continuous determination of batch average product temperature ( $T_b$  and  $T_p$ )
- Continuous determination of:
  - $R_p$ : product resistance to drying
  - $\ell$ : product dry layer thickness





# Determination of Primary and Secondary Drying Endpoints



Data spikes: TDLAS data recorded during MTM-based SMART experiment

**SMART FD Cycle**  
**5% sucrose formulation**  
**3 mL fill**  
**20 mL vials**  
**112 vials**


# Pilot Scale FD Batch Product Temperature Average



Weighted average thermocouple temperature (for edge and center vials) was calculated for comparison with TDLAS  $T_b$ . EV: Edge vials; CV: Center vials; TC: Thermocouple.  $K_v$   $2.90 \times 10^{-4}$  cal/sec.cm<sup>2</sup>.K,  $N=1620$ ,  $A_v$ : 7.17 cm<sup>2</sup>,  $\Delta H_s$ : 660 cal/sec.

$K_v$  scaled from lab FD measurements: emissivity & edge vial ratio


# Continuous Determination of Product Resistance and Ice Thickness



## TDLAS Application in QbD

- PAT tool providing key data for QbD based drying cycle development by determination:
  - Vial heat transfer coefficient ( $K_v$ )
  - Product dry layer resistance ( $R_p$ )
  - FD capability limits: Onset of choked flow
- Assessment of drying heterogeneity: prediction of # of vials completing 1<sup>o</sup> drying
- Applicable to all freeze dryer sizes enabling scale up experiments and technology transfer

# Construction of Design Space




## TDLAS Determination of:

- Kv: vial heat transfer coefficient
- Rp: product resistance to drying
- FD equipment limit: choked flow



# Kv Determination With Changing Chamber Pressure



Kuu, W.Y., Nail, S.L., Sacha, G., Rapid Determination of Vial Heat Transfer Parameters Using Tunable Diode Laser Absorption Spectroscopy (TDLAS) in Response to Step-Changes in Pressure Set-Point During Freeze-Drying, J Pharm Sci, 98(3) 2009.



# Sample Kv Data Using TDLAS Measurement

- Use of TDLAS enables generation of this table in one experiment
- Gravimetric approach requires one experiment per pressure level (12 experiments)

Pressure (mT)	$K_v$ (j/hr-cm <sup>2</sup> - °K)
25	3.58
50	5.25
75	6.11
100	7.24
125	8.20
150	9.01
175	9.75
200	10.31
250	11.21
300	12.07
350	12.92
400	13.77

Adapted from Nail & Kessler:  
Experiences with TDLAS at Laboratory & Production Scales, Garmisch 2010


## Scale up to Production Freeze Dyer Opportunity for use of TDLAS

- Lab scale QbD-based cycle development: knowledge & design space
  - Determination of vial heat transfer coefficient,  $K_v$
  - Determination of product resistance to drying,  $R_p$
  - Establishment of FD equipment limitation: choked flow measurements


# Production Scale QbD Based Cycle Modifications

- Commercial scale QbD-based cycle development: knowledge & design space
  - Evaluate freeze dryer limitation: choked flow measurements
  - Adjust lab scale dryer  $K_v$  for production dryer
    - Scale by differences in ratio of center/edge vials and shelf & wall emissivities
  - Re-evaluate design space using adjusted  $K_v$  and lab scale  $R_p$
  - Freeze dry demonstration batch using modified cycle design with data from batch ( $R_p$ )
    - Verify design space with measured values:  $dm/dt$  &  $T_p$
    - Confirm new cycle design with second demonstration batch

PC with build in high frequency source




**Quartz based sensor**, operating on the principle of temperature dependent resonance: after excitation by a modulated microwave signal (2.4 GHz) the sensor keeps on oscillating in a temperature dependent frequency. Overlaying the sensor response with the carrier signal leads to a frequency shift from which the product temperature  $T_b$  can be derived.





**1** Lab FD with TEMPRIS LAB

**2** Pilot-FD with TEMPRIS LAB




**3** Production-FD with TEMPRIS PRO










## In Summary

- Process control and PAT tools facilitate the QbD approach and conforms to regulatory expectations
- ControLyo™ technology is gaining more commercial manufacturing adoption in the industry
- Controlled nucleation minimizes variabilities in product quality
- Use of TDLAS can provide product information which can impact quality attributes



## Line of Sight across the SP Range



## References

- ICH Q8 (R2)Guidance for Industry – Pharmaceutical Development
- ICH Q11 - Development and Manufacture of Drug Substance (Chemical Entities and Biotechnological/Biological Entities)
- FDA Draft Guidance - Advancement of Emerging Technology Applications to Modernize the Pharmaceutical Manufacturing Base Guidance for Industry
- FDA Guidance - PAT — A Framework for Innovative Pharmaceutical Development, Manufacturing, and Quality Assurance
- FDA – Guide to Inspection of Lyophilization of Parenterals

## Ian Whitehall Chief Marketing Officer

